

Package ‘AssayCorrector’

May 9, 2017

Type Package

Title Detection and Correction of Spatial Bias in HTS Screens

Version 2.0.0

Date 2017-05-09

Author Bogdan Mazoure (McGill University, Department of Mathematics and Statistics), Iurie Caraus (UQAM, Departement d'Informatique), Robert Nadon (McGill University, Department of Human Genetics and Genome Quebec Innovation Center) and Vladimir Makarenkov (UQAM, Departement d'Informatique)

Maintainer Bogdan Mazoure <bogdan.mazoure@mail.mcgill.ca>

Description (1) Detects plate-specific spatial bias by identifying rows and columns of all plates of the assay affected by this bias (following the results of the Mann-Whitney U test) as well as assay-specific spatial bias by identifying well locations (i.e., well positions scanned across all plates of a given assay) affected by this bias (also following the results of the Mann-Whitney U test); (2) Allows one to correct plate-specific spatial bias using either the additive or multiplicative PMP (Partial Mean Polish) method (the most appropriate spatial bias model can be either specified by the user or determined by the program following the results of the Kolmogorov-Smirnov two-sample test) to correct the assay measurements as well as to correct assay-specific spatial bias by carrying out robust Z-scores within each plate of the assay and then traditional Z-scores across well locations.

Depends R (>= 3.3.0),
RColorBrewer,
lattice,
latticeExtra,
kSamples,
RVAideMemoire

License GPL-3

LazyData TRUE

RoxygenNote 6.0.1

R topics documented:

correct_bias	2
create_assay	2
detect_bias	3
m	4
plate7	4
plot.assay	5
print.assay	5

Index**6**

correct_bias	<i>Correct the bias present in the assay, previously detected by the detect_bias() method</i>
--------------	---

Description

correct_bias() (1) uses either the additive or multiplicative PMP (Partial Mean Polish) methods (the most appropriate spatial bias model can be either specified or determined by the program following the results of the Kolmogorov-Smirnov two-sample test) to correct the assay measurements if the plate-specific correction is specified; (2) carries out the assay-specific correction if specified.

Usage

```
correct_bias(assay, method = NULL, alpha = 0.05, type = "PA")
```

Arguments

assay	The assay to be corrected. Has to be an assay object.
method	NULL:autodetect (default), 1:additive, 2:multiplicative
alpha	Significance level threshold (defaults to 0.05)
type	P:plate-specific, A:assay-specific, PA:plate then assay-specific, AP:assay then plate-specific

Value

The corrected assay (assay object)

Examples

```
assay<-create_assay(m)
detected<-detect_bias(assay)
corrected<-correct_bias(detected,method=2)
```

create_assay	<i>Create a new assay</i>
--------------	---------------------------

Description

create_assay makes a new object of class assay. You should pass this object to detect_bias() and correct_bias() methods

Usage

```
create_assay(m, ctrl = NA)
```

Arguments

<code>m</code>	The assay you want to be corrected
<code>ctrl</code>	An optional boolean array of the same dimensions as <code>m</code> . Each entry is 1 if the well is a control well, 0 otherwise. All control wells are excluded from all computations

Value

`assay` The created assay object. It contains the following fields:

`n` The HTS matrix of raw measurements

`ctrl` The binary matrix of control wells

`biasPositions` The binary matrix where 1:well is biased, 0:well is unbiased, as suggested by Mann-Whitney test

`mCorrected` The HTS matrix of corrected measurements, initialized to a zero array, and subsequently storing the corrected version of `m` via `correct_bias()`

`biasType` Vector of length `p`, where `p` is the number of plates. It tells, for each plate of the assay, A:Additive trend, M:Multiplicative trend, U:Undetermined trend and C>Error-free plate.

`biasModel` Vector of length `p`, where `p` is the number of plates. It tells, for each plate of the assay, the most likely spatial bias model (1 through 6)

`biasConf` Vector of length `p`, where `p` is the number of plates. It tells, for each plate of the assay, the confidence in the model, (0 - lowest to 3- highest).

Examples

```
# Fictive 8x12x5 assay
assay<-create_assay(m)
# Plate 7 taken from Carralot et al. 2012
assay<-create_assay(plate7)
```

`detect_bias`

Detect the type of bias present in the assay

Description

`detect` (1) identifies rows and columns of all plates of the assay affected by spatial bias (following the results of the Mann-Whitney U test); (2) identifies well locations (i.e., well positions scanned across all plates of a given assay) affected by spatial bias (also following the results of the Mann-Whitney U test).

Usage

```
detect_bias(assay, alpha = 0.01, type = "P", test = "AD")
```

Arguments

assay	The assay to be corrected. Has to be an assay object.
alpha	Significance level threshold (defaults to 0.01)
type	P:plate-specific, A:assay-specific, PA:plate then assay-specific, AP:assay then plate-specific
test	KS:Kolmogorov-Smirnov (1933), AD:Anderson-Darling (1952), CVM:Cramer-von-Mises (1928)

Value

The corrected assay (assay object)

Examples

```
assay<-create_assay(m)
detected<-detect_bias(assay)
```

m	<i>Artificially generated HTS assay</i>
---	---

Description

A 8x12x5 assay (class array) that was generated with partial column bias in the last column. Gaussian noise was then added

Usage

```
m
```

Format

An object of class array of dimension 8 x 12 x 5.

plate7	<i>Plate 7 of raw measurements obtained by Carralot et al. 2012 (Institut Pasteur Korea)</i>
--------	--

Description

The seventh plate take from an experimental assay, which showed the evidence of existence of multiplicative kind of bias in screening technologies.

Usage

```
plate7
```

Format

An object of class matrix with 16 rows and 20 columns.

plot.assay	<i>Plot assay plate-wise</i>
------------	------------------------------

Description

plot.assay plots a hit map of the assay (only one plate at a time)

Usage

```
## S3 method for class 'assay'  
plot(x, ..., plate = 1, type = "R")
```

Arguments

x	The assay you want to plot
...	Ellipsis to be passed to the levelplotplot() function from the lattice package
plate	The plate number (Default:1)
type	Either "R" - raw assay, "C" - corrected assay (if it exists) or "P" - spatial bias position

Value

None

print.assay	<i>Print assay summary</i>
-------------	----------------------------

Description

print.assay simply prints a summary of the HTS assay

Usage

```
## S3 method for class 'assay'  
print(x, ..., plate = 1)
```

Arguments

x	The assay you want to print
...	Ellipsis to be passed to the default print() function
plate	The plate number (Default:1)

Value

None

Index

*Topic **datasets**

m, [4](#)

plate7, [4](#)

correct_bias, [2](#)

create_assay, [2](#)

detect_bias, [3](#)

m, [4](#)

plate7, [4](#)

plot.assay, [5](#)

print.assay, [5](#)